# lumie® Arabica™

**OPERATING INSTRUCTIONS**

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**BRUKSANVISNING**

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English

## OPERATING INSTRUCTIONS

### Arabica

---

- 1 Precautions
  - 2 Safety
  - 3 Getting started
  - 4 How long?
  - 5 When?
  - 6 How often?
  - 7 When should I start?
  - 8 Problem solving
  - 9 Replacement bulbs
-

English

## INTRODUCTION

Congratulations on buying Arabica. Using Arabica every day will provide you with extra bright light during the winter months and can help alleviate SAD (Seasonal Affective Disorder) and the Winter Blues. Light therapy is a safe way to get the bright light your body needs and thousands of people have benefited, feeling more energetic and with a brighter outlook.

Please read these instructions carefully to get the most out of your Arabica.

*J Affect Disord.* 2000 Jan-Mar;57(1-3):55-61

*Am J Psychiatry.* 2005 Apr;162(4):656-62

# 1 Precautions

Research shows that light therapy is very safe. There are no problems with ultraviolet (UV) light as Arabica emits only a tiny amount, far less than normal daylight. However, please check with your doctor first if any of the following apply:

- you have a major eye problem
- you have been advised to avoid bright light
- you are taking medicine where the instructions advise you to avoid bright light
- you have, or have had, deep depression\*

\*This is especially important if you are also using medication, since Arabica is designed to help improve your mood.

If you have any problems whilst using Arabica, or if you've been using the light for a week but haven't noticed any effects, please call the Lumie Careline (01954 780500) or speak to your doctor.

**IMPORTANT:** The guidelines in this leaflet are not medical directions for treating any condition. If your doctor has suggested using light therapy, follow their instructions and discuss any effects with them

## 2 Safety

- Keep the Arabica away from water and damp. For indoor use only.
- Keep the cable away from heated surfaces.
- Children using the Arabica, or playing near it, should be supervised.
- Always use the correct bulb (24W/PLL/860). Replacements are available from Lumie.

## 3 Getting started

- Choose a place where it won't slide about or get knocked over.
- Position your Arabica at arm's length (50cm/20inches) with the plastic window pointing towards your face.
- Plug in and switch on.

Although you don't need to stare at it, the light must reach your eyes to have any effect. This means you can use Arabica whilst you are working at the computer, watching TV, reading, etc. Dark glasses, tinted lenses or falling asleep will all limit the effects of light therapy. You should notice positive benefits after using Arabica for 3 or 4 days in a row.

## 4 How long?

Everybody's different but, if you use Arabica about 50cm/20inches from your face as suggested, between 1 and 2 hours (cumulative – not necessarily in a single session) is usually enough. If you're more comfortable with it further away, then less light will reach your eyes and you'll need to use it for longer. If you find your eyes ache slightly, try sitting further away or switch off for a while.

## 5 When?

It depends on your sleep cycle. If you find it difficult to wake up in the morning, switch on as soon as you get up; on the other hand, if you fall asleep too early, use Arabica in the late afternoon or early evening. If you've no time in the mornings, or can only manage a short session, then it's fine to top-up in the afternoon or early evening. Once you've found a pattern that's right for you, try to use it around the same time every day. It's best not to use the light late in the evening as you may find it difficult to get to sleep.

## 6 How often?

Once you notice the difference, you'll get into the habit of using your Arabica daily. Most people can skip a day here and there, or even the occasional long weekend, as long they resume afterwards. Your body will soon tell you if you've missed some light therapy sessions and you should use your Arabica for 3 or 4 consecutive days to make them up.

You can monitor your progress with the questionnaire in this booklet.

## 7 When should I start?

It's best to tackle the Winter Blues as soon as possible. Look out for early signs - often in September or October - and if you start to feel lethargic, anxious or irritable then start using your Arabica. If you learn to use the light early on then you should never dread the winter months again.

## 8 Problem solving

### **The light is less bright or doesn't come on at all.**

- One or both of the bulbs may have failed. Switch off and unplug the lightbox. Let it cool if necessary.
- Remove the plastic screen by gently pulling it away from the top of the unit; it unhooks at the bottom.
- Check the internal fuse – if it's blown then this will be visually evident. Replace this with a 20mm cartridge fuse type F3.15AL250V.
- Check that the bulbs are firmly in the sockets. If you have to adjust them, plug in the lightbox and switch on afterwards to see if this has worked.
- If the bulbs are still not working you will need to replace them. To remove a bulb, push the red button down hard and slide the bulb out of the socket.
- Gently slide the replacement bulb into the socket until it clicks into position.
- Do NOT remove the reflector as there are no user-serviceable parts under the reflector.

**I've replaced the bulbs but it's still not working.**

- Try replacing the internal fuse (F3.15AL250V; 5x20mm cartridge type). However, if your Arabica keeps blowing fuses, there is a problem and you should contact us.
- It is also possible that the plug fuse might have blown. The correct fuse is a 3A fuse to BS1362.

If you have any other problems with your Arabica, please call the Lumie Careline (01954 780500) before sending the unit back for repair.

## 9 Replacement bulbs

The bulbs will slowly lose their brightness so, to maintain maximum output, we recommend you change them after about 500 hours of use. Replacement bulbs are available from Lumie.

English

## Warranty

This unit is covered under warranty against any manufacturing defect for 3 years from the date of purchase. This is in addition to your statutory rights. If you are returning your Arabica for service, it must be in its original packaging including internal pieces. Please call us for further advice.

## Contact us

Lumie is a trademark of Outside In, Europe's leading specialist in light therapy. Our products are based on the latest scientific studies and we focus solely on light therapy and its benefits. If you have any suggestions that you feel could help us improve our products or service, or if there are problems you think we should know about, please call the Lumie Careline on 01954 780500 or email [info@lumie.com](mailto:info@lumie.com).

Lumie  
3 The Links, Trafalgar Way  
Bar Hill  
Cambridge  
CB23 8UD  
UK

English

## Technical specification

Subject to change without notice.

### **Lumie Arabica**

230Vac 50Hz 55W

Lamp 2x24W/PLL/860

Operating temperature 0°C-35°C

Class IIa Medical Device

---

Lumie is a trademark of Outside In (Cambridge) Limited  
incorporated in England and Wales.  
Registered Number 2647359  
VAT No GB 880 9837 71

© Lumie 2009

---

10 | [www.lumie.com](http://www.lumie.com)

## English

This questionnaire is designed to help you assess how much the lights are helping. It is based on the standard tool used by doctors. Not everyone will have every symptom but your total score should go down over the weeks. You might like to talk it over with your doctor or another person in your household.

Compared with when you are well, how have you felt during the past week?

0 = no different    1 = a little bit    2 = somewhat    3 = quite a lot    4 = definitely, badly

DEPRESSION SYMPTOMS CHART					
I have been...	before using the light	after 1 week	after 2 weeks	after 3 weeks	after 4 weeks
down and depressed					
less interested in doing things					
less interested in sex					
less interested in eating					
losing weight					
difficulty getting to sleep					
sleeping badly, waking at night					
waking too early					
feeling a failure					
feeling that life is not worth living					
tense, irritable, worried about things					
convinced I'm ill, excessively worried about my health					
thinking or speaking more slowly					
fidgety, restless, can't keep still					
feeling cut off from what's going on, as if in a dream					
paranoid, suspicious					
preoccupied, easily distracted					
physically tired e.g. heavy limbs, headaches, aching back/muscles					
other physical symptoms e.g. stomach problems, sweating, urinating frequently, dry mouth, cramp, heart palpitations, hyperventilating					
morning is worse than evening					
evening is worse than morning					
<b>Standard depression score</b>					
feeling unsociable, avoiding people					
gaining weight					
wanting to eat more than usual					
ACTUALLY eating more than usual					
craving sweet and starchy foods					
sleeping too much, feeling more tired than I should be					
worse in the afternoon or evening, then better at least an hour before bedtime					
<b>Supplemental depression score</b>					
<b>TOTAL standard and supplemental score</b>					

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Outside In (Cambridge) Ltd